


Comparison PageWriter TC Series	TC70	TC50	TC30
			
Screen Size	15 inches	10 inches	6.5 inches
Warranty	1 yr on-site	2 yr bench with loaner	3 year bench
Keyboard	Full size	¾ size (laptop)	¾ size (laptop)
Battery Life	Requires 2 (50 ECGs per battery)	2nd battery opt. (30 ECGs per battery)	2nd battery opt. (30 ECGs per battery)
16-lead PIM	Optional	Optional	Not available
Paper exit	Print from side	Print from front	Print from front
Screen movement	Closes flat	Closes flat	Tilts 18 degrees
Disclosure time	20 minutes	10 minutes	5 minutes
Event markers	15 events	10 events	6 events
Manual orders	Standard	Standard	Optional
Networked Orders (200)	Standard	Optional	Optional
Archive 200 ECGs stored	Standard	Standard	Optional
Interpretation	Standard	Standard	Optional
Disclosure	Standard	Standard	Optional
Last ECG & Query	Standard	Optional	Optional