


PageWriter TC70



Reveal more

Our DXL 16-Lead Algorithm incorporates right heart and posterior wall information across a broad range of adult and pediatric conditions.

It's as easy as 1-2-3



Just touch it


Take ECGs from the large touch screen, the keyboard, or the Patient Interface Module with one tap of the green button.

STEMI-CA

Statements suggest which coronary artery may be occluded.

Up-to-date

The nomenclature embedded into PageWriter TC70 is consistent with AHA/ACC recommendations.



Color-coded lead maps


Unique LeadCheck software tests for 19 different lead reversals to help you get it right the first time.

Critical Values

Quickly identify patients that need urgent care in support of Joint Commission Patient Safety goals.

ST Map

At a glance get a clear indication of ST elevation for quick triage.



Synchronize time


Auto set the PageWriter time with your hospital time master to obtain accurate documentation of your patient's clinical history.



Download orders

Get Orders

Import ECG orders with complete patient information from Philips TraceMasterVue.



One-touch workflow

PageWriter TC70 can be configured to automatically print, save, transfer, and retrieve a previous ECG—with one touch of a button.